



CERNER FLASH

06/03/2018

Prosolve-Finalizing Reports

In Prosolve, doctors must associate the order to the report prior to finalizing the report.

The screenshot shows a dialog box titled "Order Numbers" with a close button (X) in the top right corner. The dialog contains the following text:

This Exam performed on 2018-06-02 09:43
Please choose the appropriate
check in number from the list below

Placer Number	Effective Date	Universal Service ID
349492364	2018-06-02 08:38	MCLNM73^NM Myocardial Perfusion SPECT - Multi^
349481072	2018-06-02 08:03	Stress Nuclear Treadmill^Stress Nuclear Treadmill^

At the bottom of the dialog are two buttons: "OK" and "Cancel".