


# Cerner Flash

06/11/2023

## PathNet – Troponin Instructions

The reflex rule from the Initial Troponin order that triggered the Troponin 2 hour automatically (as the system had prior) is not yet functioning correctly in the new PathNet system. While the team works on fixing that rule, please see below for how to order Troponins to receive results as expected.

1. Initial Troponin will provide one Troponin. Please ensure a STAT collection priority for rapid results. It will not reflex the 2 hour Troponin until the rule is fixed.
2. Troponin 2 hour and 4 hour draws should be changed to have a Timed collection priority with a specific date and time or a STAT priority with an specific start date and time.
3. If the Initial Troponin and a Troponin 2 hour are ordered at the same time and the ordering provider does not adjust the start time on the 2 hour, both labs will print on the same collection label and cause staff confusion leading to a cancelled 2 hour test or two results from the same time instead of 2 hours apart.



Laboratory

Troponin T High Sens... Order 6/11/2023 12:07 EDT Plasma, 6/11/2023, Once, Nurse collect, Order for future visit

Details for **Troponin T High Sensitivity (hs-cTnT) 2 hour**

Details Order Comments Diagnoses

\*Specimen Type: Plasma

\*Requested Start Date/Time: 06/11/2023 1207 EDT

Duration:

Nurse collect:  Yes  No

\*Collection Priority:

\*Frequency:

Duration unit:

Order for future visit:

Stat

Timed Study